


Certificate of Analysis Cannabinoids

Reference:	CUR6B-152	Client:	Plantoflife
Sample date:	-----	Sample ID:	17300618
Bloomday:	-----	Sample material:	oil
Description:	Plant of Remedy		
Further information:	Turmeric Oil 30% CBD		

Abbr.	Substance	Result	unit
P-GEW	Sample weight	4.209	g
T-CBD	Total Cannabidiol (CBD + CBDA)	29.11	% (w/w)
CBD	Cannabidiol	29.11	% (w/w)
CBDA	Cannabidiolic acid	ND**	% (w/w)
T-THC	Total Tetrahydrocannabinol (THC + THCA)	ND**	% (w/w)
D9THC	D9-Tetrahydrocannabinol	ND**	% (w/w)
THCA	Tetrahydrocannabinolic acid	ND**	% (w/w)
D8THC	D8-Tetrahydrocannabinol	ND**	% (w/w)
T-CBG	Total Cannabigerol (CBG + CBGA)	ND**	% (w/w)
CBG	Cannabigerol	ND**	% (w/w)
CBGA	Cannabigerolic acid	ND**	% (w/w)
CBN	Cannabinol	ND**	% (w/w)
CBC	Cannabichromene	ND**	% (w/w)
THCV	Tetrahydrocannabivarin	ND**	% (w/w)
CBDV	Cannabidivarin	0.10	% (w/w)
CBDVA	Cannabidivarinic Acid	ND**	% (w/w)

Picture of the received sample on 10/08/2022



Head of Laboratory Services



Ing. Christian Fuczik, Chemist
Analysis reviewed - last changes:
12/08/2022 at 10:17

Footnote:

** ND =not detectable. The measured value was below the limit of detection of 0.01 % or 100 mg/kg.

The expected measurement uncertainty varies with substance and concentration and can be assumed to be a maximum of 5 %.

For the calculations of the equivalent sums, the respective acid forms were multiplied by the factor 0.877 or 0.878 to conclude the equivalent amount of the neutral form.

Method of analysis: HPLC-DAD (High Performance Liquid Chromatography - Diode Array Detector) according to Ph.Eur. 2.2.29 (European Pharmacopoeia)

This Certificate of Analysis may only be reproduced as a whole and not in parts. Any alteration is punishable under § 223 StGB (Austrian Penal Code) (forgery of documents).